




Power Column, Flatwall Headwall System and Contour Headwall System





P976-02 Infusion Support Arm

- Designed to hold two transfer poles (P158A thru Hill-Rom) and can be located almost anywhere on the track
- SWL 80 lbs



P984-00 Infusion Pump Holder

- Four hangers, requires two P970A01 universal holders



P972-01 Post-Type Utility Mount

- Requires P970A01 Universal Holder
- Length of post - 2 1/4"



P973-01 Utility/Bottle Slide Track Mount

- Requires P970A01 Universal Holder



P1020D1 Utility/Bottle Slide Flush Mount



P971H00 Slide-Mounted IV Hanger

- Outer Dimensions: 2"W x 4 1/2"H
- Inner Dimensions: 1 3/4"W x 4"H
- Requires P973-01 or P1020D1 Utility/Bottle Slide Track Mount
- Requires P970A01 Universal Holder

	Power Column	Flatwall	Contour	Wall Mount
P976-02 Infusion Support Arm	●	●		
P984-00 Infusion Pump Holder	●	●		
P972-01 Post-Type Utility Mount	●	●		
P973-01 Utility/Bottle Slide Track Mount	●	●	●	
P1020D1 Utility/Bottle Slide Flush Mount	●	●	●	
P971H00 Slide-Mounted IV Hanger	●	●	●	



P158A Transfer Pole

- Includes four IV hangers
- Retracted Position - 32"
- Fully Extended - 57 1/4"
- Diameter of top pole - 1/2"
- Diameter of outer pole - 1"
- Requires P976-02 Infusion Support Arm



P972L01 Post-Type Long Utility Mount

- Requires P970A01 Universal Holder
- Length of post - 12"



P970A01 Universal Holder

Power Column	Flatwall	Contour	Wall Mount
●	●		
●	●		
●	●	●	

**P962CMU Hewlett-Packard Merlin
CPU/Satellite Rack Mount**

- Optional mounting of CPU on back of Power Column

P006638 Monitor Mount for the side of Power Column. .

**P129244x01 Monitor Mount 3 point for
Contour**



P100702 Monitor Slide, Wall Mount Only

	Power Column	Flatwall	Contour	Wall Mount
P962CMU Hewlett-Packard Merlin CPU/Satellite Rack Mount	●			
P006638 Monitor Mount for the side of Power Column. .				
P129244x01 Monitor Mount 3 point for Contour				
P100702 Monitor Slide, Wall Mount Only				●



P17041201X01 Aneroid Sphygmomanometer with coil and a standard size cuff

- The Tycos aneroid sphygmomanometer mount, cuff basket and universal holder are ordered separately.



P975-01 Aneroid Mount and Cuff Basket for 5091-38

- Requires P970A01 Universal Holder

P975A01 Aneroid Mount for P17041201X01 WELCH ALLYN 767



P975A01 Aneroid Mount for a Welch Allyn Model 7670-01 .

- Requires a P970A01 Universal Holder



P973-S2 Integrated Diagnostic System Mount/Welch Allyn Model 767

- Requires P970A01 Universal Holder



P973-S1 Integrated Diagnostic System Mount/Welch Allyn Model 747

- Requires P970A01 Universal Holder

	Power Column	Flatwall	Contour	Wall Mount
P17041201X01 Aneroid Sphygmomanometer with coil and a standard size cuff	●	●	●	
P975-01 Aneroid Mount and Cuff Basket for 5091-38	●	●	●	
P975A01 Aneroid Mount for P17041201X01 WELCH ALLYN 767				
P975A01 Aneroid Mount for a Welch Allyn Model 7670-01	●	●	●	
P973-S2 Integrated Diagnostic System Mount/Welch Allyn Model 767	●	●	●	
P973-S1 Integrated Diagnostic System Mount/Welch Allyn Model 747	●	●	●	



P951W03 Large Pivoting Storage Basket

- Inner dimensions - 9 9/16" W x 9 9/16" H x 12 7/8" D
- Functional loading up to 24 lbs.
- Light neutral color

P17010302X01 Basket 6 1/2" W x 5" H x 4" D Requires a P970A01 Universal Holder



P951S03 Small Pivoting Storage Basket

- Inner dimensions - 9 9/16" W x 9 9/16" H x 4" D
- Functional loading up to 4 lbs.
- Light neutral color

P951S03X01 Small Basket with Dividers



P952S01 Small Storage Basket

- Inner dimensions 21 3/8" W x 4 3/8" H x 2 7/8" D
- Functional loading up to 3 lbs
- Light Neutral Color

P952S01X0001 Small Basket for 15" Power Column



P952L01 Large Storage Basket

- Inner dimensions 20 7/16" W x 8 13/16" H x 6 7/8" D
- Functional loading up to 15 lbs
- Light Neutral color basket with Black bumper



P672E03 Storage Basket for 19" Flatwall

- Inner dimensions: 18 11/16" W x 7" H x 9 1/16" D
- Functional loading up to 15 lbs



P672E04 Storage Basket for 28" Flatwall

- Inner dimensions: 27 11/16" W x 7" H x 9 1/16" D
- Functional loading up to 20 lbs
- Light Neutral Color with Black Bumper

Power Column	Flatwall	Contour	Wall Mount
●	●	●	
●	●	●	
●			
●			
	●		
		●	



P986S01BMO7 Small Swivel Shelf

- Light neutral laminate
- Overall dimensions - 15" x 12"
- Weight capacity - 50 lbs.



P986-03BMO7 Large Swivel Shelf

- Light neutral laminate
- This swivel shelf can be rotated 360°.
- Overall dimensions - 15" x 22"
- Weight capacity - 50 lbs.

P961E02X11 Adjustable Shelf for 15" Power Column



P961E02BMO7 Adjustable Shelf for Power Column (shown) .

- Light neutral laminate
- Overall dimensions - 22 1/8"W x 12 1/8"D
- Weight capacity - 50 lbs.

P673E21BMO7 Adjustable Shelf for 19" Flatwall

- Light neutral laminate
- Overall dimensions - 22 3/8"W x 12 1/8"D
- Weight capacity - 50 lbs.

P673E30BMO7 Adjustable Shelf for 28" Flatwall

- Light neutral laminate
- Overall dimensions - 29 3/8"W x 12 1/8"D
- Weight capacity - 50 lbs.

P673E16B Adjustable Shelf for Contour 15 1/2 " W x 12"D.



P674E01 19" Flatwall Guard/Bumper (shown)

Also Available:

- **P674E02 28" Flatwall Guard/Bumper**
- **P674E03 37" Flatwall Guard/Bumper**
- **P953-00 23 1/8" Power Column Guard/Bumper**
- **P953-01 15 1/8" Power Column Guard/Bumper**
- **P674E07 Contour Guard/Bumper**



P674E05 28" Flatwall Bumper with Three Bottle Racks (shown) ...

Also Available:

- **P674E04 19" Flatwall Bumper with Two Bottle Racks**
- **P674E06 37" Flatwall Bumper with Four Bottle Racks**

	Power Column	Flatwall	Contour	Wall Mount
P986S01BMO7 Small Swivel Shelf	●	●	●	
P986-03BMO7 Large Swivel Shelf	●	●	●	
P961E02X11 Adjustable Shelf for 15" Power Column				
P961E02BMO7 Adjustable Shelf for Power Column (shown)	●*			
P673E21BMO7 Adjustable Shelf for 19" Flatwall		●		
P673E30BMO7 Adjustable Shelf for 28" Flatwall		●		
P673E16B Adjustable Shelf for Contour 15 1/2 " W x 12"D.			●	
P674E01 19" Flatwall Guard/Bumper (shown)		●		
P674E02 28" Flatwall Guard/Bumper		●		
P674E03 37" Flatwall Guard/Bumper		●		
P953-00 23 1/8" Power Column Guard/Bumper	●			
P953-01 15 1/8" Power Column Guard/Bumper	●			
P674E07 Contour Guard/Bumper			●	
P674E05 28" Flatwall Bumper with Three Bottle Racks (shown)		●		
P674E04 19" Flatwall Bumper with Two Bottle Racks		●		
P674E06 37" Flatwall Bumper with Four Bottle Racks		●		

*Not available for the 15 1/8" Power Column



P963D01 Arm Lamp - Fluorescent

- Arm extends to 44"
- Requires power outlet and P100803 Mount for track mounting
- 13" Cord Length



P963M00 - Wall Mount for Arm Light



P100803 - Mount Bracket



P966A00 Utility Light

- Requires power outlet
- 24" Cord Length



P1022PH Resuscitator Bag Holder






- Requires P970A01 Universal Holder



P973-B1 Bird Blender Mount

- Requires P970A01 Universal Holder

Power Column	Flatwall	Contour	Wall Mount
●	●	●	●
	●	●	●
			●
●	●	●	
●	●	●	
●	●	●	

		Power Column	Flatwall	Contour	Wall Mount
 <p>P967B00 STAT Clock•Timer Panel Mount</p> <ul style="list-style-type: none"> • For replacement only 		●	●		
 <p>P964B00 STAT Clock•Timer Pod Mount</p> <ul style="list-style-type: none"> • Hard wiring required • Requires P970A01 Universal Holder 		●	●		
 <p>P967B01 STAT Clock•Timer Wall Mount</p> <ul style="list-style-type: none"> • Requires four gang masonry box 3 1/2" deep provided by others 		●	●		
 <p>P973-C1 Cable Organizer</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 		●	●		
 <p>P958-00 Floor Mounted Docking Shoe</p>					
 <p>P1020C1 Extended Bottle Slide</p> <ul style="list-style-type: none"> • Dimensions - 1 5/8"W x 4 5/16"H • Requires accessory bar or arm to mount 		●	●	●	



P979-00 6' Vertical Wall Tracks with Keyway at Center

Wall mount only



P985A01 23" Power Column Accessory Bar

For accessories that mount on this arm, please see Integris B accessories brochure

P985A05 15" Power Column Accessory Bar

For accessories that mount on this arm, please see Integris B accessories brochure

P985A02 19" Flatwall Accessory Bar

For accessories that mount on this arm, please see Integris B accessories brochure

P985A03 28" Flatwall Accessory Bar

For accessories that mount on this arm, please see Integris B accessories brochure

P985A04 37" Flatwall Accessory Bar

For accessories that mount on this arm, please see Integris B accessories brochure

P985A06 Contour Accessory Bar

For accessories that mount on this arm, please see Integris B accessories brochure

Power Column	Flatwall	Contour	Wall Mount
			●
●			
●			
	●		
	●		
	●		
		●	

