



Latitude™ Arm System Accessories





P997L00 Monitor Arm

Provides flexible patient monitor positioning, variable height arm for light weight monitors up to 30lbs.



P962-00 Tilt Swivel Monitor Arm

- Requires monitor manufacturer name and model number.



P100803 Mount Bracket

Arm Light Mount to mount P965-01.



P1020A1 Utility/Bottle Slide

- Mounts on P985P03/P985P04



P1020B1 Bird Blender Mount

Outer dimensions: 1-15/16" W x 4-1/2" H

Inner dimensions: 3/4" W x 2-7/8" H

- Mounts on P985P03/P985P04



P1022PH Resuscitator Bag Holder

- Requires P970A01 Universal Holder



P103601 Cable Organizers

- Mounts on P985P03/P985P04



P158A Transfer Pole

- Includes four IV hangers
- Retracted Position - 32"
- Fully Extended - 57 1/4"
- Diameter of top pole - 1/2"
- Diameter of outer pole - 1"
- Requires P976-02 Infusion Support Arm



P951S03 Small Pivoting Storage Basket

- Inner dimensions - 9 9/16"W x 9 9/16"H x 4"D
- Functional loading up to 4 lbs.
- Light neutral color



P951W03 Large Pivoting Storage Basket

- Inner dimensions - 12 1/2"W x 12 1/2"H x 12 7/8"D
- Functional loading up to 24 lbs.
- Light neutral color



P963C01 Arm Lamp with Incandescent Bulb

Arm extends to 46". Requires power outlet and Universal Light Mount (P100803). 13" Cord length.



P963C02 Arm Lamp with Fluorescent and Incandescent Bulbs

Arm extends to 46". Requires power outlet and Universal Light Mount (P100803). 13" Cord length.



P964B00 STAT Clock•Timer Pod Mount

- Hard wiring required
- Requires P970A01 Universal Holder



P966A00 Utility Light

- Requires power outlet
- 24" Cord Length



P970A01 Universal Holder



P971H00 Slide-Mounted IV Hanger

- Outer Dimensions: 2"W x 4 1/2"H
- Inner Dimensions: 1 3/4"W x 4"H
- Requires P973-01 or P1020D1 Utility/Bottle Slide Track Mount
- Requires P970A01 Universal Holder



P1020D1 Utility/Bottle Slide Flush Mount



P972-01 Post-Type Short Utility Mount

- Requires P970A01 Universal Holder
- Length of post - 2 1/4"



P972L01 Post-Type Long Utility Mount

- Requires P970A01 Universal Holder
- Length of post - 12"



P973-01 Utility/Bottle Slide Track Mount

- Requires P970A01 Universal Holder



P973-B1 Bird Blender Mount

- Requires P970A01 Universal Holder



P973-C1 Cable Organizer

- Requires P970A01 Universal Holder



P973-S1 Integrated Diagnostic System Mount/ Welch Allyn Model 747

- Requires P970A01 Universal Holder



P973-S2 Integrated Diagnostic System Mount/ Welch Allyn Model 767

- Requires P970A01 Universal Holder



P975-01 Aneroid Mount and Cuff Basket

- Requires P970A01 Universal Holder




P976-02 Infusion Support Arm

- Designed to hold two transfer poles (P158A) and can be located almost anywhere on the track



P978-01 Aneroid Sphygmomanometer with coil and a standard size cuff

- The Tycos aneroid sphygmomanometer mount, cuff basket and universal holder are ordered separately.



P984-00 Infusion Pump Holder

- Four hangers, requires two P970A01 universal holders



P985P03 Accessory Arm (left-shown)

P985P04 Accessory Arm (right)

- Right Hand 18"
- Left Hand 18"

For accessories that mount on these arms, please see Horizon/Uniflex Accessories brochure



P986-03BMO7 Large Swivel Shelf

- Light neutral laminate
- This swivel shelf can be rotated 360°.
- Overall dimensions - 15" x 22"
- Weight capacity - 50 lbs.



P986S01BMO7 Small Swivel Shelf

- Light neutral laminate
- Overall dimensions - 15" x 12"
- Weight capacity - 50 lbs.



P17000501 Short Post

2-1/4" Post, 3/4" Diameter Includes Universal Holder (P17000601). Functional loading up to 15 lbs.

- Mounts onto P662A0701



P17000601 Universal Holder

Designed to accept various accessories.

- Mounts onto P662A0701



P17000701 Universal Pole Holder

(C-Clamp) Accepts poles from 3/4" to 1" in diameter. Designed to accommodate the P158A ISS Transfer Pole.


- Mounts onto P662A0701



P17010302 Basket

6-1/2"W x 4"D x 5"H Functional loading to 5 lbs. (Basket mounts in front of rail.)

- Mounts onto P662A0701



P17010402 Basket

11"W x 4"D x 4"H Functional loading to 5 lbs.
(Basket mounts in front of rail.)

- Mounts onto P662A0701



P17010501 Basket - Pivoting

9-9/16"W x 9-9/16"D x 4"H Functional loading to 5 lbs.
Designed to pivot. Inner liner is provided to contain supplies.

- Mounts onto P662A0701



P17043001 Conchatherm Humidifier & Heater Mount

Slide Adapter only.
• Mounts onto P662A0701



P17043101 Puritan-Bennett Humidifier Holder

For Cascade humidifier. Post is 4" high, 1-1/2" diameter.

- Mounts onto P662A0701



P17070101 Standard Slide Adapter

1-3/16" inner width.
• Mounts onto P662A0701



P17070201 Bird™ Blender Slide

Used to mount oxygen blenders, bottle holders, suction canisters and other slide-mounted devices. Adapter for Bird Blenders.

- Mounts onto P662A0701



P17071001 Dual Clip Assembly

- Mounts onto P662A0701



P17072001 Cable Organizer

- Mounts onto P662A0701



P17073001 Vertical Cord Wrap

- Mounts onto P662A0701



P17073101 Horizontal Cord Wrap

- Mounts onto P662A0701



P662A0701 Horizontal Accessory Rail

WittRock
Latitude™ Arm System
Accessories



WittRock reserves the right to make changes without notice in design, specifications and models. The only warranty WittRock makes is the express written warranty extended on the sale or rental of its products.

© 2005 WittRock
ALL RIGHTS RESERVED.

140189 Rev. 1 09/12/05